

[W&M ScholarWorks](https://scholarworks.wm.edu/)

[VIMS Articles](https://scholarworks.wm.edu/vimsarticles) [Virginia Institute of Marine Science](https://scholarworks.wm.edu/vims)

1988

Determination Of Tributyltin In The Marine-Environment

RG Huggett Virginia Institute of Marine Science

Michael A. Unger Virginia Institute of Marine Science

FA Espourteille Virginia Institute of Marine Science

CD Rice Virginia Institute of Marine Science

Follow this and additional works at: [https://scholarworks.wm.edu/vimsarticles](https://scholarworks.wm.edu/vimsarticles?utm_source=scholarworks.wm.edu%2Fvimsarticles%2F1491&utm_medium=PDF&utm_campaign=PDFCoverPages)

Part of the [Environmental Chemistry Commons](http://network.bepress.com/hgg/discipline/134?utm_source=scholarworks.wm.edu%2Fvimsarticles%2F1491&utm_medium=PDF&utm_campaign=PDFCoverPages)

Recommended Citation

Huggett, RG; Unger, Michael A.; Espourteille, FA; and Rice, CD, "Determination Of Tributyltin In The Marine-Environment" (1988). VIMS Articles. 1491. [https://scholarworks.wm.edu/vimsarticles/1491](https://scholarworks.wm.edu/vimsarticles/1491?utm_source=scholarworks.wm.edu%2Fvimsarticles%2F1491&utm_medium=PDF&utm_campaign=PDFCoverPages)

This Article is brought to you for free and open access by the Virginia Institute of Marine Science at W&M ScholarWorks. It has been accepted for inclusion in VIMS Articles by an authorized administrator of W&M ScholarWorks. For more information, please contact [scholarworks@wm.edu.](mailto:scholarworks@wm.edu)

Accuracy in Trace Analysis

III. Quantitation in Material, Environmental, Clinical, and Nutrient Analysis

Environmental Analysis

Determination of Tributyltin in the Marine Environment

R. G. Huggett, M. A. Unger, F. A. Espourteille, and C. D. Rice

Virginia Institute of Marine Science School of Marine Science College of William and Mary Gloucester Pt., VA 23062

Introduction

Tributyltin (TBT) is a biocide used in antifouling paints to protect hulls from nuisance organisms such as barnacles, worms and algae. The use of TBT paints has increased over the past decade due to its effectiveness as an antifoulant which is related to its toxicity. Water concentrations of less than 100 ng L^{-1} have been shown to harm some aquatic species in laboratory tests and observations in the natural environment indicate that levels below 10 ng L^{-1} may be harmful. Tributyltin is bioconcentrated by many species to levels of one thousand, or more, times ambient water concentrations. Sediment-water partitioning coefficients for TBT of 100-10,000 have been reported [1]. The extreme toxicity of TBT challenges the analytical chemist to accurately and precisely determine ambient TBT concentrations in water at or below I ng L^{-1} and in sediments and tissue at concentrations ranging from μ g kg⁻¹ to mg kg⁻¹.

Methods and Procedures

Tributyltin is often derivatized to convert the substance to a more volatile compound to facilitate analyses [2,3]. Our methodology for water involves liquid-liquid extraction of an acidified sample with tropolone and n-hexane. This is followed by the generation of n-hexyltributyltin, using the Grignard reagent, n-hexyl magnesium bromide. The

derivatized extract is then cleaned by fluorisil column chromatography and analyzed by capillary gas chromatography with flame photometric detection or gas chromatography-mass spectrometry. Tripentyltin chloride is added as an internal standard before extraction of the sample [4].

The methodology for quantification of TBT in sediments and tissue differs from the water procedure only in the extraction steps. Anhydrous sodium sulfate and precipitated silica are added to wet tissue or sediment to dessicate the sample. It is then frozen to lyse cells and ground to a fine powder consistency. The TBT is removed from the sample by soxhlet extraction with n-hexane. The extract is then processed by the same procedure used for water extracts [5].

Discussion

The methodologies described above are straightforward, but a unique set of problems emerge when attempting to quantitate TBT at low ng L^{-1} concentrations. For instance, most commercial supplies of the Grignard reagent contain organotins as contaminants. This occurs because the same inert atmosphere apparatus is often used to synthesize organotins and Grignard reagents. Some commercial Grignard reagents were found to contain TBT at levels as high as several mg kg^{-1} . A contamination-free commercial source of n-hexyl magnesium bromide has since been found which now allows us to quantify TBT in water at less than 1 ng L^{-1} .

Sediment and tissue analyses are more complex than those for water. The material is difficult to extract from these matrices and the analyst must be cautious not to degrade or alter the substance in the process. Additionally, compounds that may interfere with TBT analyses are more abundant in tissue and sediment than in water. This is particularly true for samples collected from marinas and harbors which experience hydrocarbon pollution from oil spills and bilge waters.

Accuracy in Trace Analysis

Sampling frequency is an important aspect of any analytical investigation involving chemical contamination of the natural environment. Our work with TBT confirms this statement. The natural variability of TBT concentrations in water, sediment or tissue is large. Table 1 shows the concentrations of TBT in 19 individual oysters, *Crassostrea virginica,* collected from a single location in the James River, VA [6]. The lowest concentration differs from the highest by more than a factor of two. Water samples collected at the same location and at the same phase of the tide can vary by an order of magnitude within a one week period. Samples collected at the same time but a few hundred meters apart can be vastly different in concentration [7]. There are a number of factors responsible for this natural variability. These include wind conditions which influence both the amplitude of the tides and the mixing of the water mass, the ever-changing number of vessels in the vicinity of the sampling sites and the age of the paint films, since the release of TBT from the paints is somewhat time dependent. Figure I shows the concentration of TBT in surface water (30 cm below surface) collected from the same location on high slack tides over an 18-month period. Such data indicate that TBT water monitoring programs, and perhaps monitoring programs for

Table 1. Concentrations of tributyltin in oysters, *Crassostrea virginica,* from a single location in the James River, VA

Oyster	Tributyltin concentration ^a μ g kg ⁻¹
1	301
	285
$\frac{2}{3}$	211
4	515
5	280
6	345
7	386
8	254
9	374
10	258
11	411
12	446
13	480
14	311
15	264
16	379
17	373
18	325
19	325
	$x \pm sd = 343 \pm 81$

a Concentrations given as TBT, dry weight.

other dissolved constituents, must be carefully designed with particular attention given to natural variability if short-term trend analyses (e.g., over I or 2 years) are desired.

The proper interpretation of environmental chemical monitoring data depends not only on the accuracy and precision of the analyses but also on how well the samples represent the system being investigated. The authors suggest that more attention should be given to sampling in the total analytical scheme.

Acknowledgments

The authors wish to acknowledge the laboratory support of D. J. Westbrook and E. J. Travelstead. Financial assistance was provided by the Virginia State Water Control Board, the U.S. Navy and the Virginia Institute of Marine Science.

Figure 1. Tributyltin in near surface waters from Sarah Creek, VA over an 18-month period.

References

- [1] Unger, M. A., MacIntyre, W. G., and Huggett, R. J., Equilibrium Sorption of Tributyltin Chloride by Chesapeake Bay Sediments, In: Proceedings of the Oceans '87, International Organotin Symposium, Halifax, Canada (1987).
- [2] Maguire, R. J., Chau, *Y.,* Bengert, G. A., Hale, E. J., wong, P. T. S., and Kramer, O., Environ. Sci. Tech. 16, 698 (1982).
- [3] Matthias, C. L., Olson, G. J., Brinckman, F. E., and Bellama, J. M., A Comprehensive Method for the Determination of Aquatic Butyltin and Butylmethyltin Species at Ultra-Trace Levels Using Simultaneous Hybridization/Extraction with G. C./F.P.D., presented in part at American Chemical Society, Miami Beach, FL (1985).
- [4] Unger, M. A., MacIntyre, W. G., Greaves, J., and Huggett, R. J., Chemosphere 15, 461 (1986).

Accuracy in Trace Analysis

- [5] Rice, C. D., Espourteille, F. A., and Huggett, R. J., Appl. Organometallic Chem. 6, in press.
- [6] Espourteille, F. A., Tributyltin in Oysters, *Crassostrea virginica,* and Sediments from the Chesapeake Bay, Masters Thesis, School of Marine Science, College of William and Mary, Gloucester Pt., VA (in preparation).
- [7] Huggett, R. J., Unger, M. A., and Westbrook, D. J., Organotin Concentrations in the Southern Chesapeake Bay, In: Proceedings of Ocean '86 Organotin Symposium, Marine Technology Society, Washington, DC, 1262-1265 (1986).

Analysis of Gaseous and Particle-Associated PAH and Nitroarenes in Ambient Air

Janet Arey, Barbara Zielinska, Roger Atkinson, and Arthur M. Winer

Statewide Air Pollution Research Center University of California Riverside, CA 92521

The organic species present in ambient air vary considerably in volatility and concentration and are distributed between the gaseous and particulate phases. This distribution depends on, among other factors, ambient temperature and the composition of the particulate phase. The major sources of anthropogenic chemicals emitted into the atmosphere are combustion and volatilization, and the more volatile species are often present in high concentrations. This is demonstrated in figure I for a series of aromatic hydrocarbons. Literature values for ambient benzene and toluene concentrations [1] and the ambient concentration data we obtained for polycyclic aromatic hydrocarbons (PAH) at Glendora, CA, in August 1986 are shown plotted against their room temperature vapor pressures [2,3]. Thus, to measure ambient concentrations of the complete range of PAH and nitroarenes (which are generally at least an order of magnitude less abundant than the parent PAH), complementary sampling and analysis techniques are required.

Indicated in table I are the sampling technique(s) which provided quantitative concentration valuesfor the listed PAH during 12-h sampling periods. The sampling techniques employed were: sampling at 1 L min⁻¹ onto cartridges containing 0.1 g of Tenax GC solid adsorbent (low-flow Tenax), sampling at 10 L min⁻¹ onto 0.6 g of Tenax (high-flow Tenax), sampling at 30 SCFM onto four polyurethane foam (PUF) plugs located in series downstream from Teflon-impregnated glass fiber (TIGF) filters in modified high-volume samplers, and standard high-volume sampling at 40 SCFM onto TIGF filters. Crucial to the quantification

Table 1. PAH quantitatively measured and the sampling techniques employed at Glendora, CA, in August 1986

m.w.	PAH	Low- flow Tenax	High- flow Tenax	PUF plugs	Filter	Comment
128	Naphthalene	x				Ave. 2% on back-up Tenax
142	1-Methylnaphthalene X		х			Low- and high-flow Tenax agree
142	2-Methylnaphthalene X		X			Low- and high-flow Tenax agree
154	Biphenyl		х			Too low to quantify on low-flow Tenax
154	Acenaphthene		x			Too low to quantify on low-flow Tenax
166	Fluorene		X			Σ PUFs ~25% high-flow value
178	Phenanthrene		X	X		$<$ 10% on 4th PUF
178	Anthracene			x		Too low to quantify on high-flow Tenax
184	Dibenzothiophene			x		
202	Fluoranthene			X		Partially gas-phase
						Partially "blown-off" filter
202	Pyrene			X		$<$ 10% on filter
						Partially gas-phase
						Partially "blown-off" filter
						$< 10\%$ on filter
228	Chrysene/Triphenylene				ΣΧ	20–50% present on PUFs
	Higher Molecular Weight PAH				x	